

NAD ASSAY

REAGENTS

- 1: GLYCYL-Glycine BUFFER (x 250 ml) at pH 7.4
 - 0.065 M glycyl-glycine in distilled water
 - 3.05 g Niacinamide
 - 6.05 ml 95% ethanol

- 2: ALCOHOL DEHYDROGENASE (ADH)
 - from a stock solution of 29.4 mg/ml pick up 20 μ l and dilute in 1 ml of phosphate buffer (10 mM at pH 7). **20 μ l/ml**

- 3: MYF
 - 2.5 mg/ml distilled water

- 4: PHENAZINE METOSULFATE (PMS)
 - 4 mg/ml distilled water

METHOD

1: SAMPLES PREPARATION

In 12 well plates:

- Remove supernatant off cells
- Add 250 μ l Perchloric acid 0.5 N and scrape the cells
- Transfer the cells to eppendorf tubes
- Add 70 μ l of KOH 3M to neutralize
- Centrifuge for 2 minutes at 40 °C

2: ASSAY

- In 96 well plates add for each well:
 - 190 μ l Glycyl-glycine buffer
 - 10 μ l ADH
 - 10 μ l MTT
 - 20 μ L PMS
- Incubate the plate for 10 minutes at 37C in the plate reader
- Add 20 μ l of sample to each wells with the above reagents
- Make the blank with all the above reagents but without the samples
- Make the standard with all the above reagents with 20 μ l of NAD 1 μ M
- Read at 560 nm at 0, 10 and 20 minutes at 37 °C.

3: CALCULATION

- Subtract blank from the OD values of samples (or read subtracting the blank)
- Plot OD values on y axis and time on x axis and calculate the simple slope
- Divide the obtained slope by the standard slope
- Data are expressed as nmol/ml or n cells of NAD

CELL CULTURE PREPARATION FOR ATP (6-well and 12-well plate prep.)

CELL SCRAPING

1. After removal of the cell culture plates from the incubator, keep on ice until each well is ready for scraping.
 2. Aspirate 500 μ l media from each well for nitrate and place in a well of a 96-well plate. (store at -20 C until nitrate analysis)
- NOTE: omit step #2 if nitrate determination is not needed.
3. Aspirate remaining media from wells without touching cell monolayer.
 - 4a. For 6-well plates: Immediately add 1.0cc of iced perchloric acid (0.6N) to each well.
 - 4b. For 12-well plates: Immediately add 0.5cc of iced perchloric acid (0.6N) to each well.
 5. With a cell scraper, scrape cells from each well, gently mixing with the perchloric acid.
 6. Transfer cells into labelled eppendorf tubes and place on ice.
 7. Vortex each epp. tube for 10sec. at speed 5.

CELL SONICATION

8. Using Virtis membrane disruptor, sonicate each epp. tube containing the cell extract for 10 sec at setting #5. (position the epp. tube in an ice bath while sonicating)
9. Leave epp. tubes on ice for 60min.

NEUTRALIZATION

- 10a. Vortex epp tube 5sec
- 10bi. For 6-well plates: Immediately add 450 μ l of "ICED" K₂HPO₄ (1.0M pH 12)
- 10b2. For 12-well plates: Immediately add 225 μ l of "ICED" K₂HPO₄ (1.0M pH 12)
- 10c. Vortex again for 5sec
11. Centrifuge at 10,000g for 10 min at 4 C

FILTRATION

- (Using 3cc syringes, 19ga X 1.5inch needles, gelman acrodisc or equiv. 0.2 μ syringe filters.)
12. Aspirate supernate from the pellet with the 3cc syringe and needle.
 13. Remove needle and attach syringe filter to the syringe.
 14. Filter sample into a labelled micro-vial.
 15. Immediately drop vial into container of liquid nitrogen to freeze instantly.
 16. Store at -70 C until ATP analysis

DIGESTION OF CELL PELLETT

17. Invert the pellet to allow to dry at room temp.
18. Add 1.0cc of NaOH (1.0N) to pellet and vortex until pellet is mixed.
19. Store at room temp. until protein is assayed.

PROTEIN (Pellet) ANALYSIS

1. Centrifuge the protein digest for 10 minute at 4,000RPM at room temp.
2. Set-up 1 2X75mm test tubes for the protein determination.
3. In each tube, add 800 μ l dist. H₂O
4. Add 200 μ l Bio-Rad dye reagent and vortex.
5. Next add protein digest sample: add 50 μ l.
6. Immediately VORTEX again right after adding the protein sample!!
7. Transfer 200 μ l of each test tube sample into 96 well plate.

8. Prepare SpectraMax for protein analysis:

- a) Double click MacHD
- b) Click on Open Settings
- c) Click on Settings Protein/Bio-Rad
- d) Go to Template
- e) Paint the plate according to where your sample wells are located.
- f) Click on Group and go to Unknowns and Click O.K.
- h) Open drawer on plate reader SpectraMax machine and put 96 well plate in position.
- i) Go back to computer and click on READ!!!
- j) Go to file and Print Plate 1, print 2 copies for Alvin.